
MATTRESS



15 & 17 cm неіднт



Description

CORE High Resilience polyurethane foam 40 kg/m³, cut into waffle-style studs on the "sleeper" side for the prevention of bedsores, multi-articulated on the other side

Available in HR55 $kg/m^{\scriptscriptstyle 3}$ for overweight patients – bariatric bed

FABRIC M1 polyure thane coated polyester mesh cover, breathability 400 g/m² per 24 hours, green, 150 g/m².

FINISH Removable cover with a waterproof reverse zip down one length and one width on bottom edge with identification label on the "sleep" side. No need to cover the zip to ensure perfect hygiene and easy cleaning.

MAINTENANCE Machine wash at 90°C, tumble dry up to 60°C, use so apy water, very diluted bleach. Do not use ether or fabric softener

TREATMENT antibacterial, anti-dust mite and anti-fungal treated fabric.

FIRE STANDARD Permanent MI fire rated fabric. HR foam M4 fire rated

COMPLIANCE Mattress compliant with GPEM/CP 2005 (standard NF EN597 parts 1 and 2 - recommendation D1 90 GPEM/CP) and with article U23 of January 2005 Class B

COMFORT INDEX Resistance to compression 40% (kPa): 3 (40 kg/m 3), or 3.7 (55 kg/mt3)

CUSTOMER BASE Mattress for patients with a Norton score \geq 14 with low risk or stage 1 bedsores.

GUARANTEED against premature sagging.

OPTION Bevelled corners





